



# Q2 Special Offers

April - June 2026

Buy any 4 Protect+ or Eco+ Sterilisation Pouches, get another box FREE\*!

## Protect+ Self-Seal Sterilisation Pouches

Easy to open pouches with triple seal end for extra protection

- All individual Protect+ Self-Seal Sterilisation Pouches contain a sign and date area for an easy auditing process
- Manufactured from 60g medical weight sterile packaging paper and transparent blue PET/PP-plastic composite laminate film
- Both paper and plastic are fully recyclable, refer to local/national recycling regulations



### Packs of 200

Use	Size	Code
Burs & Endo Instruments	57 x 127mm	PERF-0025000
Burs & Endo Holders	85 x 165mm	PERF-0025001
Hand Instruments & HP's	85 x 253mm	PERF-0025002
Multiple Hand Instruments	135 x 278mm	PERF-0025003
Push Bar Instrument Tray	190 x 330mm	PERF-0025004
Full Instrument Tray	305 x 420mm	PERF-0025006
Full Instrument Tray	300 x 450mm	PERF-0025009

## Eco+ Sterilisation Pouches

- Colour indication to show sterilisation process has taken place
- Medical grade paper
- Bulk pack for economy

### Box 500

Size	Code
89 x 254mm	PERF-0025030
140 x 280mm	PERF-0025031
190 x 355mm	PERF-0025032



\*Mix & Match, cheapest item free

**Buy any ImpressPLUS Putty, Wash, Clear or Bite Registration material, get Pack 25 Impression Trays with Handle FREE**



## ImpressPLUS Putty

- Dark Blue High Viscosity VPS Putty
- Excellent Dimensional Stability
- High tear strength and elasticity
- 300ml base, 300ml catalyst and two colour coordinated scoops

Set Type	Working Time	Time in mouth	Code
Regular	2 min	2.30 min	PERF-0101300
Fast	1 min	2 min	PERF-0101302



## ImpressPLUS Wash

- VPS material with excellent tear strength
- Light body yellow, medium body buff and monophas purple
- 2 x 50ml cartridges, 6 disposable tips and 6 intra oral tips

Regular Set	Fast Set
Working Time: 2.00 min	Working Time: 1 min
Time in the mouth: 2.30 min	Time in the mouth: 2 min

Set Type	Code Regular Set	Code Fast Set
Light Body	PERF-0101360	PERF-0101362
Medium Body	PERF-0101370	PERF-0101372
Monophas Body	PERF-0101380	PERF-0101382



## ImpressPLUS Bite Registration

Thixotropic light blue VPS material for bite registration, a quick die for temporary inlays and onlays, copying transfers and primary impressions with bite trays

- 2 x 50ml cartridges and 6 disposable tips

Set Type	Working Time	Time in mouth	Code
Regular	1.30 min	3 min	PERF-0101350
Fast	0.40 min	1 min	PERF-0101352
Super Fast	0.30 min	0.45 min	PERF-0101354



## ImpressPLUS Clear

Transparent VPS A-Silicone material

- Ideal for creating transparent matrices for highly aesthetic restorations
- Can be used alongside injectable composites to create highly aesthetic restorations such as indirect veneers
- 2 x 50ml cartridges & 6 tips



Set Type	Working Time	Time in mouth	Code
Fast	1 min	2 min	PERF-0100348

## Plastic Impression Trays with Handles

Perforated Impression Trays with moulded plastic handles

### Pack 25

Uses	Code
Upper Large	PERF-0113030
Lower Large	PERF-0113031
Upper Medium	PERF-0113032
Lower Medium	PERF-0113033

Uses	Code
Upper Small	PERF-0113034
Lower Small	PERF-0113035
Partial UR/LL	PERF-0113036
Partial UL/LR	PERF-0113037
Lower Anterior	PERF-0113038





## SealEDPLUS

### Pit & Fissure Sealant

- A light cure radiopaque fluoride releasing sealant designed to protect the fissures.
- An opaque low viscosity and flowable material
- The fluoride releases during polymerisation and for several days after application providing maximum protection against caries
- Fluoride release of 2.52ppm after 24 hours

Product	Code
2 x 2g Syringe & 5 Tips	PERF-0101215



## BridgelT

### Temporary Crown & Bridge Material

- A self-curing material offering excellent aesthetic result combined with optimal mechanical properties making it ideal for provisional restorations
- The formulation ensures low heat generation during the intra oral setting protecting against damaging the pulp
- High flexural strength allows for easy placement and the compressive strength provides exceptional strength and resistance of the final restoration against masticatory forces

#### Technical Data

Mixing ratio	1:1
Working time	60 seconds
Curing time in mouth	1.00-2.00 mins
Flexural strength (MPa)	95
Compressive strength (MPa)	298

#### 1:1 Ratio Cartridge System 75g and 6 Mixing Tips

Shade	Code
A2	PERF-0101610
A3	PERF-0101611



## RestorePLUS LC Cavity Liner

### Light curing, glass reinforced cavity liner

- Provides insulating layer to reduce risk of ongoing sensitivity
- Glass reinforcement enhances mechanical strength, wear resistance and durability
- Fluoride releasing to prevent secondary caries - 1.80ppm after 24 hours
- High compressive and flexural strength
- If applying to exposed pulp, place PerfectaLINE Calcium Hydroxide Liner over pulp initially

Product	Code
2 x 2g Syringes & Tips	PERF-0100027



## Protect+ Instrument Cleaning Sachets

- Individual pre-dosed, soluble cleaning sachets, specially formulated with surface-active cleaner
- Removes dirt, dried blood and other organic matter
- Pre-dosed, plastic-free, fully soluble sachets
- 1 x 17.5g sachet makes up to 3L of solution
- Unperfumed and suitable for use in Ultrasonic units

Product	Code
Box of 50 x 17.5g Sachets	PERF-0020541

## ErgoPLUS

Premium range of Deluxe 9.5mm round hollow handled instruments designed for superior ergonomics and comfort.

The range of Scalers and curettes are manufactured with premium "stay sharp" material for optimum strength and sharpness.

The lightweight handle helps minimise user fatigue, reducing potential strain and injury.



### Probes

#### Single Ended

9.5mm Handles approx. 18g weight

No.	Code
6	PERF-0152010
14W	PERF-0152016



#### Perio Screening Probes

9.5mm Handles approx. 18g weight

No.	Code
CPITN-C	PERF-0152015
UNC15	PERF-0152017



#### Double Ended

9.5mm Handles approx. 18g weight



No.	Code
DG16 Endo Explorer	PERF-0152040

Was	Now	Was	Now
€17.00	€12.75	£14.62	£10.95

### Mirror Handle




Code	Was	Now	Was	Now
PERF-0152000	€8.20	€6.15	£7.07	£5.30

### Scalers

#### Double Ended Hollow

9.5mm Handles approx. 18g weight



No.	Code	Was	Now	Was	Now
H6/H7	PERF-0152120	€23.99	€17.99	£20.63	£15.45
204S	PERF-0152121				

### Gracey Curettes

#### Double Ended Hollow

9.5mm Handles approx. 18g weight



No.	Code	Was	Now	Was	Now
1/2	PERF-0152100	€23.00	€17.25	£19.80	£14.85
3/4	PERF-0152101				
5/6	PERF-0152102				
7/8	PERF-0152103				
9/10	PERF-0152104				
11/12	PERF-0152105				
13/14	PERF-0152106				

Ask your local Sales Consultant for an ErgoPLUS Demonstration

Follow DMI on Social Media



@DMILtd



Dental-Medical-Ireland



@dmi\_dental

Scan the QR to see full range of Perfection Plus products



Republic of Ireland



Northern Ireland

**Terms & Conditions:** Prices and discounts are subject to change without notice and are ex-warehouse. Value Added Tax will be added at the prevailing rate. Insurance, handling, packing and delivery charges are amounts which the customer shall pay in addition. Orders are accepted with the understanding that Goods will be billed at prices prevailing on date of shipment. Payment terms are 30 days from date of invoice for consumables. Goods approved for return by the Company, which are received in a saleable condition, must be returned within 30 days of invoice date, otherwise a handling charge of 20% plus vat will apply. Goods returned to the Company are done so at the Customer's risk and must be shipped prepaid. A goods returned docket form and a copy of the original invoice must accompany the return. Offers valid until 30th June 2026, whilst stocks last.

**Contact Details:** Call: Dublin: 01 427 3700/Cork: 021 602 0544/Lisburn: 028 9260 1000. Email: info@dmi.ie. Web: www.dmi.ie/www.dmi.co.uk