



DENTAL MEDICAL IRELAND

# Dental Stools

## 2023

Relax knowing we have the  
right chair for you



The ergonomic features on your stool can profoundly influence your musculoskeletal health for years to come. The ideal stool should support your back, neck, and shoulder posture.

DMI offer a wide range of doctors and nurses' stools from leading manufacturers which are designed to flex, conform and adjust to your every move.

The dental stool should be the most important chair in your practice. The amount of time a dentist/dental assistant spends in the chair can vary. Therefore, a good dental stool must provide comfort as well as durability and be simple to function.



EXCLUDES VAT E&OE

Call: 028 9260 1000 Web: [www.dmi.co.uk](http://www.dmi.co.uk) Email: [info@dmico.uk](mailto:info@dmico.uk)

Unless stated otherwise, all saddle chairs have a two-part seat, Ø 560 mm Salli Base and standard Ø 65 mm castors. Available also with Ø 540 mm base. Height adjustment by gas spring. The upholstery is genuine leather, artificial leather, ESD artificial leather, fabric. Maximum load 120 kg.

### 1. Salli SwingFit

- Swing mechanism makes it possible for the user to exercise while sitting
- The width of the seat can be adjusted
- The tightness of the swing mechanism can be adjusted
- 6 cm narrower Small available



Standard: Salli Base



Small

640 £

### 2. Salli MultiAdjuster

- Seat width and tilt adjustments
- 6 cm narrower Small available



Ø 540 mm base



Small

602 £

### 3. Salli Swing

- Swing mechanism makes it possible for the user to exercise while sitting
- The tightness of the swing mechanism can be adjusted



545 £

### 4. Salli Twin

- Seat with tilt adjustment



508 £

### 5. Salli Classic

- Undivided seat with tilt adjustment



521 £

### 6. Salli Light / 7. Salli Slim

- Single black polyurethane piece with a groove in the middle
- Salli Slim is for children, and for women who want a smaller seat
- Maximum load 100 kg
- Plastic base (Ø 580 mm for Light, Ø 480 mm for Slim) as a standard\*
- Finnish design



LISWTCHA  
LITCHA  
LINTCHA



SLSWCHA  
SLNTCHA  
SLTCHA

Basic  
202 £

Tilt  
227 £

Swing  
267 £

### 8. Salli Premium

- Salli SwingFit with your own name
- Black Salli Base
- Black Premium castors
- Not available as Small model
- Not available in ESD



674 £

\*Price with aluminium base (Ø 540 mm) +29 £  
Salli Base +50 £





Salli offers ESD (= electrostatic discharge) protected alternatives for electronics industry and office work. ESD chairs come with an antistatic upholstery, a zinc electroplated seat frame, and Ø 65 mm ESD castors. With Swing, SwingFit (picture), Twin, MultiAdjuster, Classic. Accessories: Elbow Rest, Stretching Support, Top & Bottom.

## 9. Salli Strong

- Extra strong structure, bigger seat
- Maximum load 150 kg
- Tilt and height adjustment in the same lever
- Base Ø 600 mm
- Short, medium or long special gas spring, metallic
- No accessories can be attached



620 £

SALLISTRO

## 10. Salli Expert

- Back and elbow supports
- Height adjustment by hand and foot
- Base Ø 540 mm



Salli Twin **1968 £**  
Salli MultiAdjuster **1058 £**  
Salli Small MultiAdjuster **1058 £**

EXTWTCHAI  
EXMUTCHAI  
EXSMALLMUTCHAI

Package price for Swing/SwingFit  
(when ordered afterwards)

Base Ø 600 mm included **281 £**

## 11. Allround

- Support for arms in front or on the sides, also behind the chair for stretching the back
- With SwingFit, Swing, MultiAdjuster, Twin or Classic
- Attached on the gas spring, curve length approx. 61 cm, width approx. 8.5 cm
- Upholstery leather, artificial leather or fabric
- Ø 600 mm base with Swing/SwingFit needed



250 £

ALLAROUND

## 12. Stretching Support

- Stretching Support is meant for momentary stretching only, not as a backrest
- With Classic, MultiAdjuster, Twin, Light Tilt or Slim Tilt
- Upholstery leather, artificial leather, ESD artificial leather, fabric or black polyurethane
- Ø 600 mm base needed



108 £

STRSUP

\*Leather (picture), polyurethane upholstery **-37 £**

ACCESSORIES

**salli**  
SIT HAPPY



Salli small multiadjuster with an ergorest support:


- The seat provides adjustable gap and tilt
- Gives dynamic support for the arms



# BAMBACH SADDLE SEAT



- Improves Sitting Posture
- Relieves Pain
- Reduces Fatigue
- Reduces Stiffness In Joints
- Increases Mobility
- Improves Circulation



## Colours Available



### The Bambach Saddle Seat without Back

- Seat Tilt
- Height adjustable gas stem
- Black 5 star base (chrome base optional)
- Standard furniture twin wheel hooded castors
- Available in a variety of colours
- Ask for details

### The Bambach Saddle Seat with Back

- Seat Tilt
- Height adjustable gas stem
- Black 5 star base (chrome base optional)
- Standard furniture twin wheel hooded castors
- Available in a variety of colours
- Ask for details



**Try before you buy  
14 day free trial**

Call for details

**Bambach - £535**

**Bambach with Back rest - £680**

- Prevents and cures spinal strain
- Provides and keeps the natural sitting position of pelvis and spinal column
- Adjustable in height and tiltable seat surface, pelvis keeps its upright position, also when bending forward
- Seat surface is covered with genuine English leather
- 12 months guarantee
- Available in a variety of colours

# Tronwind Stools



## Colours Available



MF1-Burgandy	MF2-Purple	MF3-Red	MF4-Orange	MF5-Apricot
MF6-Blue	MF7-Sky	MF8-Cyan	MF9-Green	MF10-Light Green
MF11-Brown	MF12-Curry	MF13-Light Brown	MF14-Wheat	MF15-Creamy White
MF16-Black	MF17-Charcoal Grey	MF18-Grey	MF19-Grey White	MF20-White

## Tronwind Nurse's Stool

- Arm swing support tube 360 degrees around the seat
- Ratchet mechanism allowing armrest to move closer and further away from the operative
- Oil-less bearing connect between arm swing support tube and gas lift
- Hand stitch micro fibre leather
- High density foam which ensures the long lasting shape
- High antimicrobial resistance
- Superior abrasion resistance

Code	Description	Price
TRON-TID01-MF1-6	Orange Micro Fibre Leather Nurse's Stool	£340
TRON-TD01-MF6	Blue Micro Fibre Leather Nurse's Stool	£340
TRON-TD01-MF9	Green Micro Fibre Leather Nurse's Stool	£340
TRON-TD01-MF16	Black Micro Fibre Leather Nurse's Stool	£340
TRON-TD01-MF18	Grey Micro Fibre Leather Nurse's Stool	£340



## Tronwind Doctor's Stool

- Stepless two way adjustment, allows vertical alignment for proper lumbar support and dynamic horizontal adjustment
- Hand stitch micro fibre leather
- High density foam which ensures the long lasting shape
- High antimicrobial resistance
- Superior abrasion resistance

Code	Description	Price
TRON-TD02-MF17	Charcoal Micro Fibre Leather Doctor's Stool	£340
TRON-TD02-MF6	Blue Micro Fibre Leather Doctor's Stool	£340
TRON-TD02-MF7	Light Blue Micro Fibre Leather Doctor's Stool	£340
TRON-TD02-MF16	Black Micro Fibre Leather Doctor's Stool	£340
TRON-TD02-MF18	Grey Micro Fibre Leather Doctor's Stool	£340





## Murray Gemini Stools



**MURR-SGEM/GT**



**MURR-HGEM/SAB**

- Five star aluminium base with twin wheel castors
- Chassis and metalwork finished in RAL 9002 cream
- All Murray stools are available in wide range of colours to match any dental chair
- £35 surcharge per stool to use existing chair manufacturers material to match colour.

Code	Description	Price
MURR-SGEM/GT	Dentist's stool with lumbar support back	£520
MURR-HGEM/SAB	DSA Nurse stool with foot ring and arm	£649
MURRAY GEMINI PACKAGE	Dentist's stool and Nurse's stool	£1089

## Murray Advance Stools



**MURR-SADV/GT**



**MURR-HADV/GN**

- Recently redesigned to offer very good value for money
- Tilt action and adjustable height
- Available in a variety of colours
- £35 surcharge per stool to use existing chair manufacturers material to match colour.

Code	Description	Price
MURR-SADV/GT	Dentists Stool Standard height, with tilt 430 - 560mm	£535
MURR-HADV/GN	Nurses Stool , without tilt, 420 - 710mm	£540
MURRAY ADVANCE PACKAGE	Dentist's stool and Nurse's stool	£1025

## Murray Luna Stools



**MURR-S-LUN-GT**



**MURR-H-LUN-GN**

- A quality stool which is reasonably priced, the H-LUN-GN features an adjustable high cylinder and footring. The S-LUN-GT is a Dentist stool with a double curvature seat, adjustable tilt and backrest with a standard height cylinder. £35 surcharge per stool to use existing chair manufacturers material to match colour.

Code	Description	Price
MURR-S-LUN-GT	Doctor Stool	£345
MURR-HLUN/GN	Nurses Stool	£345

## Murray Saddle Stool



**MURR-SADD-S**

- Saddle seats help to maintain an upright posture when seated. Feet are placed flat on the floor at a greater width than possible with conventional seating. The weight of the legs is taken through the feet, the upper body is stabilised. £35 surcharge per stool to use existing chair manufacturers material to match colour

Code	Description	Price
MURR-SADD-S	Standard height 450 - 630mm from floor	£515

# Dentsply Sirona Stools

Adapted to you from head to toe. Hugo, Carl or Theo – Which one will you choose?

- Ergonomic and intuitive sitting position
- Innovative design and intelligent shape for that extra bit of comfort
- Choice of two different gas spring lengths for different heights

## The color schemes

**The elegant world**  
Elegant upholstery colors for a classic, timeless look.



**The natural world**  
Soft, warm shades for a soothing treatment atmosphere.



**The vital world**  
Vivid colors and bright tones for a modern, lively atmosphere.



\* For Theo: Upholstery colors also available as Lounge upholstery.



Code	Description	Price
SIR-HUGOFREEHAND	Hugo freehand with foot release	£1,460
SIR-HUGOMANUAL	Hugo manual with manual release	£1,370
SIR-HUGOMANUALPLUS	Hugo manual plus with manual release and foot ring	£1,525
SIR-CARL	Carl manual with 360° manual release	£935
SIR-CARLPLUS	Carl manual plus with 360° manual release and foot ring	£1,025
SIR-THEO	Theo working stool	£935
SIR-PAUL	Paul manual	£559
SIR-PAULPLUS	Paul manual plus with foot ring	£645

## Gately Stools

### Dentist SRS06



The SRS 06 has a large backrest area to hug the lower back. A moulded seat for comfort. 3 lever mechanism gives height adjustment and free float motion to the back and seat if required. Each can be used independently as free float.

**£435**

GMS-SRS06N

### Nurse SRA09



This stool has a comfortable circular seating area. Adjustable seat height. The off-set padded backrest/arm/torso support give a comfortable working position. Complete with foot rest & taller gas lift.

**£399**

GMS-SRA09N

### Saddle GMS10



The saddle seating gives a seating position allowing the spine to maintain a natural 'S' shape. The saddle seating position enables you to turn and swivel while feet remain in contact with the floor maintaining in a central position which

**£485**

GMS-GMS10N

# Adec Stools



## A-dec 500 Stools

### Dynamic Support

Our feature-rich A-dec 500 stools have a dynamic seating system with unique performance zones that mold to the shape of the body. These performance zones work in unison to keep you supported while you shift, lean and work. Other features, such as the contoured backrest, cradle the spine for ultimate comfort. With A-dec's ergonomic seating solutions, you can focus on what's most important—your patients.

#### STOOL FEATURES:

##### Tilt Seat

Distributes weight from seat to feet allowing proper circulation.

##### Contoured Backrest

Cradles the spine, while encouraging blood flow. Contoured backrest is standard on the A-dec 521 doctor's stool and optional on the A-dec 522 assistant's stool.

##### Upholstery

Minimal seams make clean up easy.

##### Standard Foot Ring (A-dec 522 Assistant's stool)

Helps practitioner maintain correct posture.



**Whisper Quiet Casters**  
The design provides smooth, ultra-quiet movement on both car-



**Adjustment Controls**  
Numbered paddles are easy to access for quick adjustments and easy



**Backrest Height Control**  
A push button makes access and one handed adjustments easy.

A-dec 521 Doctor's Stool

A-dec 522 Assistant's Stool



## A-dec 400 Stools

### Iconic Durability

You deserve comfort from seating that's also going to last. For years, our durable seating solutions have been helping dental practitioners sit in a healthy supported position. Features such as the anatomically shaped backrest, adjustable torso support, rolled seat edges and ergonomic tilt, reduce disc pressure and muscle pain.

Select sewn or formed upholstery in an extensive color palette that's sure to complement your office. Make your next seating investment with A-dec—we're here to support you.

A-dec 421  
Doctor's Stool

A-dec 422  
Assistant's Stool



	A-dec 421 Doctor's Stool	A-dec 422 Assistant's Stool
Contoured backrest	•	
Contoured seat	•	•
Ergonomic tilt	•	
Adjustable torso support		•
Easy-to-reach adjustment levers	•	•
Height-adjustable foot ring		•

Note: A-dec 400 stools do not have the dynamic seating system that is standard in the 500 stool line.

### A-dec 421 Doctor's Stool Find the Right Height

To choose the right stool height, measure from hips to the feet.  
Above 32" - Choose tall cylinder  
Under 32" - Choose short cylinder



#### Sewn Upholstery (available for A-dec 500 and A-dec 400 Stools)



#### Formed Upholstery (available for A-dec 400 Stools)



Actual colors may vary. See your authorized A-dec Dealer for color samples and the most current product information.

Code	Description	Price
421	Dentist	£695
422	Assistant	£795
521	Dentist	£985
521	Assistant Armrest	£325
522	Assistant	£1,080
522	Assistant (inc backrest)	£1,365
522	Assistant Foot Pedal	£199