

Performer^{LR+}



Performer^{LR+}

FEATURES

Auto light

Program your light from your touchpad. Positions 0 and 3 turns off and positions 1 and 2 turn on.



Deluxe Touchpad

Upgrade to the deluxe touchpad for added ancillary integrations.



Left or right handed, you decide

Fully ambidextrous without the use of tools, in seconds.



Assistants instrumentation

Central vacuum canister allows for easy left right positioning.



Hidden Utilities

All internal services and utilities hidden away in a sleek modern floorbox.



Built to last.

With Performer **LR**⁺, you're assured comfort and simplicity. It's part of what makes A-dec A-dec.




Designed to give you flexibility, plus access and comfort. Performer **LR**⁺ answers the need for less complexity, lower cost and exceptional choices for every pound you invest.



Upgrade to the halogen 3-axis A-dec 6300 Light. The 6300 dental light provides a balanced light pattern with ideal contrast for a clear view of tooth topography. A neutral daylight colour temperature allows accurate colour/shade matching and tissue analysis.




The delivery system can be located in either left handed or right handed positions in a matter of seconds. It also features a 300 touchpad, which offers more functionality and the option to upgrade to the deluxe version.




Legendary A-dec reliability starts with our design philosophy:- create innovative, feature-rich products without complication. The result? Easy-to-use equipment that lasts longer and requires less service.


Upholstery choices




ash




steel




black




taupe




tawny




saffron




sapphire




chianti




paprika




wedgewood




indigo




sky blue



dusty jade



sage



lemongrass

Add a splash of colour?

Paprika or Lemongrass or

Sky Blue or Chianti?

Why choose A-dec

With almost 50 years experience of producing dental equipment, A-dec works with customers in more than 100 countries ranging from single-handed and group practices through to dental schools and hospitals. Here are just some of the reasons why so many professionals choose A-dec.

Quality and reliability

A-dec has earned a reputation over many years for creating probably the most durable and reliable equipment available. It is a reputation due largely to our dedication to quality as well as clever yet simple engineering design, which also leads to ease of equipment maintenance.

Ergonomics and efficiency

A-dec offers a very wide choice of equipment so no matter how you and your team want to work we have solutions to suit. The equipment is also designed to have outstanding flexibility of movement, so everything can be positioned to maximise ergonomics and efficiency.

Functionality and expandability

A-dec equipment can accommodate and seamlessly integrate many auxiliary items, controlled via touchpad and/or foot control. To keep your unit at the cutting edge of dentistry, the equipment design supports future upgrades and expansion.

Standard Specification Performer **LR+**

Smooth seamless upholstery

Fully ambidextrous without the use of tools, in seconds

Twin armrests

Double articulating headrest

Traditional delivery system with 300 standard touchpad

SYR, 3 x Vinyl MW, and standard Tray Holder

Performer Cuspidor

3 Position Vacuum Arm SYR/SE/HVE (Dry)

Performer Light



RRP £11,495*

*Terms and Conditions apply †All prices excludes VAT

Performer **LR+** options

6300 light option	£670.00
Satelec Scaler	£850.00
Cavitron Scaler	£1575.00
Fibre Optics (Two outlets)	£321.00
Electric Motor (Endo capable)	£2750.00

Finance Solutions

The easy way to manage your practice finances.

A clean, simple, fixed rate lease or hire purchase depending on your tax relief preference, arranged within a couple of days and sometimes within hours. It couldn't be simpler.

For more information on the Performer **LR+** call our experienced territory managers



Standard Package from:
£199 PER MONTH*



Charlie Cope
Territory Manager
Northern England, Scotland
and Ireland
Phone: 07718 786 022
Email: Charlie.Cope@a-dec.co.uk



Mark Harris
Territory Manager
Midlands, the North West and
North Wales
Phone: 07718 786 023
Email: Mark.Harris@a-dec.co.uk

Deluxe Package from:
£269 PER MONTH*



Eugene O'Malley
Territory Manager
Southern England and South Wales
Phone: 07718 786 024
Email: Eugene.O'Malley@a-dec.co.uk



Christopher Cox
Public Sector
Territory Manager
Phone: 07718 786 025
Email: Chris.Cox@a-dec.co.uk

Experience the Performer **LR+** at one of A-dec's Premier Dental Showrooms, to book your appointment please call 0800 233 285. Showrooms located in Manchester - Salford Quays, The Midlands - Nuneaton and London - Egham