

Handpiece Brochure



January - March
2024

SPECIAL OFFERS



NSK

KAVO
Dental Excellence



MK-dent®
GERMANY

For more information please call 028 9260 1000 or visit www.dmi.co.uk

Offers cannot be combined with any other promotions or discount structures. Prices exclude VAT. E&OE

TURBINES

T1 Boost Turbines

Perfection that fits in your hand

- Premium class LED turbines with ergonomic lightweight titanium sleeves
- Offering high power and optimum control at all times, the T1 range are extremely quiet during operation
- KaVo & Sirona fitting
- 3 year warranty



**SAVE
47%**

SIR-6689405
T1 Boost
(KaVo Multiflex coupling)
RRP £1,549
Only £826

T2 Boost Turbines

Quality you can rely on

- You can rely on Dentsply Sirona T2 class turbines for any situation
- The continuous high performance, the ergonomic shape together with a surface finish for strong grip, these LED turbines made of stainless steel are guaranteed to satisfy
- KaVo, Sirona, W&H & NSK fitting
- 2 year warranty



**SAVE
47%**

SIR-6689637
T2 Boost
(KaVo Multiflex coupling)
RRP £1,253
Only £668

T3 Boost Turbines

Powerful and economical

- Non optic
- Made of stainless steel, these non-fibre optic turbines have 26 watts of power and ceramic ball bearings to ensure efficient working in the long run
- Powerful lasting efficiency
- KaVo & Sirona fitting
- 1 year warranty



**SAVE
47%**

SIR-6689827
T3 Boost
(KaVo Multiflex coupling)
RRP £922
Only £492

MASTERtorque Lux M9000L

KaVo's Premium Turbine

- Powerful - 23 watts
- Quiet - Only 57dB, vibration free
- Small head for improved access and patient comfort - Height 13.1mm
- Features Direct Stop Technology (DST™)
- Four port spray
- PLASMAtec™ coating for perfect grip
- Ceramic ball bearings
- Powerful illumination. Glass rod light conductor
- KaVo & Sirona fitting



KAVO-M9000L
RRP £1,645
**£987 each or
Buy Triple Pack
@ only £746 each**

EXPERTtorque E680L

- Quiet, vibration-free operation (65 dB)

- 18 Watt power
- KaVo pushbutton chuck system (up to 30 N retention power)
- Compact head casing (height: 13.6 mm/Ø:12.5 mm)
- 4 port spray for optimal cooling
- Replaceable micro water filter
- KaVo MULTiflex coupling system
- Can be sterilised in an autoclave up to 135 °C
- KaVo & Sirona fitting
- 24 month warranty



KAVO-E680L-DUOPK KAVO-E680L
RRP £1,234
**£691 each or
Buy Triple Pack
@ only £559 each**

SMARTtorque S619L

- Quiet, vibration-free operation (65dB)
- 16 Watt power at only 2.8 bar
- KaVo pushbutton chuck system (30 N retention power)
- Compact head casing (height: 14.5 mm/Ø: 12.5 mm)
- 4 port spray for optimal cooling
- Can be sterilised in an autoclave up to 135 °C
- LUX glass rod fibre-optic
- KaVo fitting
- 12 month warranty



KAVO-S619L
RRP £793
**£486 each or
Buy Triple Pack
@ only £465 each**

Ti-Max Z900KL

- Optic
- Standard Head
- Light Weight Titanium
- 26W Power
- 4-Port Spray
- 3 Years Warranty
- Offer Fits KaVo Only



NSK-1113001
RRP £1,134
**Buy 1 get 1 FREE
@ only £567 each
while stocks last**

S-Max M900

- Optic
- Standard Head
- Stainless Steel
- 26W Power
- 4-Port Spray
- 2 Years Warranty
- Fits NSK, KaVo, Sirona, W&H, Bien Air Couplings



Fitting	Product	Product Code
Kavo	M900KL	NSK-P1258001
NSK	M900L	NSK-P1254001
Sirona	M900SL	NSK-P1262001
W&H	M900WL	NSK-P1266001
Bien Air	M900BL	NSK-P1264001

**Special Price £395
Buy 4 @ £580 each
Buy 4 @ £580 each
Buy 4 @ £580 each
Buy 4 @ £580 each**

Pana-Max 2

- Non optic
- Standard Head
- Stainless Steel
- 1 Port Spray
- Fits KaVo & NSK
- 1 year warranty



Fitting	Product Code
Kavo	NSK-P1212001
NSK	NSK-P1211001

RRP £452
**Buy 1 get 1 FREE
@ only £226 each
while stocks last**

Synea Fusion TG-98LM

- 4 port spray
- Head size Ø 11,5mm
- 19 W, 360,000 rpm
- FG burs Ø 1.6mm for bur lengths up to 25mm
- Compact glass rod, ceramic ball bearings
- Coupling system Multiflex, NSK & Rotoquick



WH-30003000

RRP £1,123
£505

Prime Line Fibre Optic HP21KL

- Standard Head, Quattro Spray
- KaVo MULTiflex Connection
- Fibre-Optic
- 25W Power
- Scratch-Resistant Titanium Coating
- German Engineered
- 2 Years Warranty




MKDE-HP21KL

**£518 each or
Buy Triple Pack
@ only £440 each**

Classic-Line Optic HC21KL

- Standard head, Quattro spray
- KaVo MULTiflex Connection
- Fibre-Optic
- 23W Power
- German Engineered
- 18 Months Warranty



MKDE-HC21KL

**£478 each or
Buy Triple Pack
@ only £406 each**

Eco-Line Fibre Optic HE21KL

- Standard head, Quattro spray
- KaVo MULTiflex Connection
- Fibre-Optic
- 24W Power
- German Engineered
- 18 Months Warranty



MKDE-HE21KL

**£391 each or
Buy Triple Pack
@ only £332 each**

CONTRA ANGLES 1:1 RATIO (BLUE)

T1 C40L



- Longer service life thanks to Diamond like Carbon coated ball bearings
- Optimal illumination for an excellent view
- 1:1 ratio has single spray
- 3 year warranty

SIR-3328437 Model: C40L
1:1 Ratio (Blue)

RRP £957
Now £510

47% OFF



T2 A40L



- Robust ceramic ball bearings
- Pleasant and easy to grip thanks to its titanium sleeves
- Optimal illumination for great visibility during treatments
- 1:1 ratio has single spray
- 2 year warranty

SIR-6376185 Model: A40L
1:1 Ratio (Blue)

RRP £800
Now £426

47% OFF



T3 E40



- The T3 is a sturdy, economical contra angle which will give many years of high performance
- Non optic
 - Sturdy, titanium sleeve
 - Ceramic ball bearings
 - Internal water
 - 1 year warranty

SIR-6392687 Model: E40
1:1 Ratio (Blue)

RRP £585
Now £312

47% OFF



MASTERmatic M20L & L68



- Transmission 1:1
- Speed range up to 40,000 rpm
- Optic
- Internal air/water coolant
- Works on electric or air motors that have a universal "E" type connection
- Plasmatec Anti-Slip Finish
- Improved angulation
- Flexibility of use when paired with different INTRA LUX Heads
- 2 year warranty

KAVO-M20L

RRP £1,248
£788 each or
Buy Triple Pack
@ only £599 each



EXPERTmatic E20L



- Quiet
- Transmission 1:1
- Cellular optic
- For latch type burs
- Speed range up to 40,000 rpm (electric motor) and 20,000 rpm (air motor)
- Internal air/water coolant
- Push button auto chuck
- 1-port spray
- EXPERTgrip™
- Universal "E" type connection
- 24 month warranty

KAVO-E20L-DUOPK
KAVO-E20L

RRP £995
£610 each or
Buy Triple Pack
@ only £478 each



SMARTmatic S20



- Non optic, no water spray
- Ceramic sliders within the clamping system for enduring and secure tool holding
- Rust-free stainless steel for long-lasting durability
- Lightweight for tireless working
- Optimum balance for precise treatment
- 1:1 Ratio
- 1 year warranty

KAVO-S20

RRP £245
£177 each or
Buy Triple Pack
@ only £170 each



Ti-Max Z25L



- Optic
- Ratio: 1:1
- 1-Port Spray
- Light Weight Titanium
- 3 Year Warranty

NSK-C1038001

RRP £1,199
£959 each
Buy 3 @ only £639 each



S-Max M25L



- Optic
- 1:1
- 1 Port Spray
- Stainless Steel
- 2 Year Warranty

NSK-C1024001

RRP £791
£633 each
Buy 4 @ only £474 each



FX 25M



- Non Optic
- 1:1
- No Spray
- Stainless Steel
- 1 year warranty

NSK-C1052001

RRP £250
Buy 1 get 1 FREE
@ only £125 each



Synea Vision WG-56LT



- Optic
- 1:1
- 1 Port Spray
- Stainless Steel
- 2 Year Warranty

RRP £1,123
£449

WH-30009000



Alegria WE-56



- Non Optic
- 1:1
- No Spray
- Stainless Steel
- 1 year warranty

RRP £345
£125

WH-30116000



CA 1:1L



- Optic
- 1:1
- 3-Port Spray
- Stainless Steel
- 2 Years Warranty

BIEN-1600384-001

RRP £611
£323



SPEED INCREASING CONTRA ANGLES 1:5 (RED)

T1 C200L



- Longer service life thanks to Diamond like Carbon coated ball bearings
- Optimal illumination for an excellent view
- 1:5 ratio has four hole spray
- 3 year warranty



SIR-3328452
Model: C200L 1:5 Ratio (Red)

RRP £1,540
Now £821

47% OFF

T2 A200L



- Robust ceramic ball bearings
- Pleasant and easy to grip thanks to its titanium sleeves
- Optimal illumination for great visibility during treatments
- 1:5 ratio has four hole spray
- 2 year warranty



SIR-6376177
Model: A200L 1:5 Ratio (Red)
SIR-6576065
AS200L (Short) 1:5 Ratio (Red)

RRP £1,079
Now £575

47% OFF

T3 E200



- Contra-Angle Red
- Ratio 1:5
- Non optic
- Spray 4N
- Ceramic bearings
- Head Diameter 8.7
- Head Height 14.6
- RPM 200,000
- 1 year warranty



SIR-6392679
RRP £866
Now £462

47% OFF

MASTERmatic M25L

- Speed increasing handpiece ideal for crown and veneer preparation
- 1:5 speed increasing ratio
- Extremely quiet - Only 55dB
- Features Direct Stop Technology
- No suck back. Replaceable spray microfilter
- Three port spray
- Plasmatec coating for perfect grip
- Ceramic ball bearings
- Powerful illumination
- 2 year warranty



KAVO-M25L

RRP £1,858
£1,065 each or
Buy Triple Pack
@ only £857 each

EXPERTmatic E25L

- Speed-increasing handpiece transmission 1:5 (red)
- Fast and precise working thanks to even smaller head and ceramic bearings
- Best reliability thanks to internally conducted spray water and air, 3-port spray, and replaceable water filter
- Perfect lighting thanks to the 25,000 Lux glass rod light conductor
- Better Grip
- 24 month warranty



KAVO-E25L-DUOPK
KAVO-E25L

RRP £1,365
£819 each or
Buy Triple Pack
@ only £655 each

Ti-Max Z95L Switch

- Optic contra-angle
- 1:5 Increasing
- Quattro spray
- Ceramic ball-bearings
- World first switch between water and spray on handpiece
- 3 year Warranty



NSK-C1034001SW

RRP £1,476
£1,181 each or
Buy 3 @ only £787 each

S-Max M95L

- Optic
- 1:5 Speed Increasing
- 4-Port Spray
- Stainless Steel
- 2 year warranty



NSK-C1023001

RRP £976
£781 each or
Buy 4 @ only £586 each

S-Max M95

- Non Optic
- 1:5 Increasing
- Max Speed 200,000 min-1
- Head Size Ø9.9 x H 14.5 mm
- For FG Burs (Ø1.6)
- Body Material Stainless Steel
- Water Spray Quattro Spray
- 2 year warranty



NSK-C1026001

RRP £791
£633 each or
Buy 4 @ only £474 each

EVO.15 1:5 L



RRP £1,316
£696

- Standard or Micro Series
- Small head (10% smaller)
- CoolTouch+™ heat-arresting technology
- Shockproof head
- 55 dBA
- Multi-strand optical glass technology
- Up to 3 years warranty

BIEN-1600940-001 Micro Series
BIEN-1600941-001 Standard

CA 1:5 L



RRP £1,097
£580

- Standard or Micro Series
- Multi-strand optical glass technology
- 58 dBA
- Cool Touch™
- Up to 3 Years Warranty

BIEN-1600690-001 Micro Series
BIEN-1600386-001 Standard

AIR MOTORS

AM1006 with LED

- Chrome coating
- LED Light
- 5,500 K colour temperature
- 20,000 rpm non-load speed
- Left and right direction of rotation
- Internal spray
- 6-pin tubing
- 1 year warranty

MKDE-AM1006

£339



181L with LED

- Speed range up to 20,000 rpm
- Forward / reverse speed adjustment
- Internal air / water coolant
- 6-pin connection
- Fits all attachments optic or non-optic with Universal "E" type connections
- 2 year warranty

KAVO-181L

£998



AM-25RM

- Non optic
- No spray
- New integrated 360° rotation
- Mid west connection
- 1 year warranty

WH-30108000

£229



AM0004

- Max. 20 000rpm
- No spray
- Standard INTRAmatic® connection
- Suitable for all common instruments
- Non fibre optic
- 1 year warranty

MKDE-AM0004

£175



STRAIGHT HANDPIECES 1:1 RATIO

Alegra HE-43

- For handpiece and contra-angle burs with Ø 2.35 mm (ISO 1797-1)
- Drive speed: 40,000 min-1
- Outer diameter of the sheath: 20 mm
- Non optic
- No spray
- 1 year warranty

WH-30119000

£120



T4 Model BH40

- Stainless steel
- Bur fastening via push button
- For use with burs with a 2.35mm diameter stem
- Max. speed: 40,000rpm
- No spray
- No optic
- 1-year warranty

SIR-6401843

£140



Smartmatic S10

- 1:1 (blue handpiece)
- Main application: Restorative
- Max. drive operating speed: 40.000
- Non optic
- No spray
- 1 year warranty

Kavo-S10

£183



FX 65M



- Non Optic
- 1:1 Straight
- No spray
- Stainless Steel
- 1 year warranty

NSK-H1014001

£250 each

Buy 1 get 1 FREE

@ only £125 each



COUPLINGS FOR KAVO TYPE TURBINES

KaVo Multiflex 465LED

- LED Optic with water control
- With 4-hole standard connections
- Rotary ring spray regulation with multiple indent locking positions
- Non-return valve (water)
- 1 year warranty

Kavo-465LED

£329



KaVo Multiflex 465LRN

- Optic with water control
- With 4-hole standard connections
- New lamp holder
- Rotary ring spray regulation with multiple indentlocking positions
- Non-return valve (water)
- 1 year warranty

Kavo-465LRN

£280



KaVo Multiflex 465RN

- With 4-hole standard connections
- Rotary ring spray regulation with multiple indent locking positions
- Non-return valve (water)
- 1 year warranty

Kavo-465RN

£180



KaVo Multiflex 460 LED

- LED Optic
- No Spray Adjustment
- For Kavo Fitting Handpieces
- 1 Year Warranty

Kavo-460LED

£182



KaVo Multiflex 460E

- Non optic
- No water control
- 1 year warranty

Kavo-460E

£122



MK-DENT Multiflex

- 6 hole quick connector with LED
- With water control
- 1 year warranty



MKDE-QC6016KW

£179



MK-DENT Multiflex

- 6 hole quick connector with LED
- No water control
- 1 year warranty

MKDE-QC6016K

£149



COUPLINGS FOR W&H TYPE TURBINES

RQ-34



- Optic (Bulb in Handpiece)
- With Spray Adjustment Ring
- For W&H Fitting Handpieces
- 1 year warranty

WH-10403400

£176



RQ-24



- Optic (Bulb in Handpiece)
- No Spray Adjustment
- For W&H Fitting Handpieces
- 1 year warranty

WH-10402400

£159



RQ-14



- Non optic
- Water control
- Free 360° instrument rotation
- Continuous regulation of the spray water supply
- Roto Quick non-retraction valve
- 1 year warranty

WH-10401400

£135



RQ-04



- Non Optic
- No water control
- 1 year warranty

WH-10400400

£101



HANDPIECE MAINTENANCE/AIR POLISHERS

Sirona LED coupling

- Optic
- 1 year warranty

SIR-6319706



Dentsply
Sirona

£324

Coupling for NSK - PTL CL

- LED optic
- No Spray Adjustment
- For NSK Fitting Handpieces
- 1 year warranty

NSK-P1001600



NSK

£186

QUATTROcare Plus 2124A

The cost-efficient solution for optimum care and long-term maintenance of handpieces.

- Cleans and oils any 4 handpieces in under 1 minute
- Only one-minute care cycle
- Program-controlled, automated care process
- Optimum lubrication, thanks to automatic, outlet-specific, metered dosing
- Easy handling
- No leakage of excess residual oil from handpiece
- Modern and ergonomic design

KAVO-1-008-3805

KAVO
Dental Excellence



RRP £2,432 **£1,702**

PROPHYflex 4 Air Polisher Box

- Easy to use, two handpieces in one – perfect for supragingival and subgingival applications - Whether for prophylaxis or periodontology, fissure sealing, aesthetic dentistry, prosthetics or orthodontics, it can be used on any teeth, any patient and almost any indication
- PROPHYflex 4 is available in three fresh and exciting colors, each available as a MULTiflex coupling or Sirona coupling option

Each box contains:

Prophy Box PROPHYflex 4 contains:

KaVo PROPHYflex 4 Wave

KaVo PROPHYpearls Neutral (80 Sticks)

KaVo PROPHYpearls Selection (80 Sticks)

Kerr OptiView Small Kit (2 pcs.)

Kavo-1-013-5078 Blue

Kavo-1-013-5079 Lime

Kavo-1-013-5080 Flamingo

KAVO
Dental Excellence

£1,172



Perio Box

£283

Enables sub gingival cleaning

Contains: 1 perio kit, 1 perio powder (4 x 100g),
1 Kerr perio probe (2pcs) Kavo-1-013-5081

RONDOflex plus 360

KAVO
Dental Excellence

- Highly efficient, clean work, light-weight and accurate
- 360 degree rotation of the tip for optimum access and visibility
- Gentle treatment without drilling, without injection, without noise, without vibrations.
- Pleasant for your patient, particularly for sensitive patients and children
- Incomparable flexibility through problem-free use at various workplaces thanks to MULTiflex coupling: simply push on and start.
- Includes 2 cannulae KAVO-01-002-2179



RRP £1,919 **£1,371**

Perio - Mate

- Used in the treatment of periodontal pockets and peri-implantitis
- The slim nozzle tip easily inserted into tight pockets and the root surfaces of adjacent teeth without causing the patient discomfort
- Powder Control Ring makes fine-tuning powder flow simple can easily be rotated using the fingertips even while performing complex operations
- Easily dismantled for cleaning and decontamination

NSK-Y1002654



RRP £1,215
£974

**Plus FREE 2nd
handpiece worth
£427**

Bundle Offer Buy Prophy - Mate & Perio - Mate together **Only £1,750**

iCare

NSK

- Cleans 4 Handpieces
- Internal cleaning and lubrication
- Chuck cleaning function ensures that the chuck is thoroughly cleaned and lubricated
- Purging select mode ensures complete removal of excess oil
- Integrated large oil tank (1.2 litre) reduces burden of refill
- 28cm (W) X 24cm (D) X 35cm (H)
- Air Line Required

NSK-Y1002787

£1,301



Prophy - Mate Neo

- Air Powered Tooth Polishing System
- Twin nozzle
- Comfort design
- Twin joints rotate smoothly 360° degrees for easy operability
- One touch operation
- Selectable 60° and 80° nozzles
- The handpiece can be washed in the thermosinfector
- Autoclavable up to 135°C
- Fits KaVo, NSK, Sirona, W&H, Bien Air

NSK



**Plus FREE 2nd
handpiece worth
£427**

RRP £1,041
£834

Universal Oil Spray

- High-quality universal care for all KaVo contra-angles, handpieces, turbines, heads and air motors
- Ensures a long life for all KaVo instruments
- Environmentally friendly propellant
- 500ml

KAVO
Dental Excellence



KAVO-0-411-9640

£28.79

Quattrocare Plus Oil Spray

- Designed to clean, lubricate and purge your instruments
- Maintenance spray for KaVo turbines, air motors, straight and contra-angle handpieces
- Only for use in the QUATTROcare handpiece oiler
- 500ml

KAVO
Dental Excellence



KAVO-1-005-4525

£34.98

Oil Nozzle Intra (Slow handpieces)

- Handpiece lubricating oil spray nozzle for slow handpieces

KAVO
Dental Excellence



KAVO-0-411-9911

£8.82

Oil Nozzle Multiflex (Fast handpieces)

- Handpiece lubricating oil spray nozzle for fast handpieces

KAVO
Dental Excellence



KAVO-0-411-9921

£8.82

Oil Spray

- For use with all types of handpiece
- Non-toxic and non-allergenic
- No odour or taste
- Good wetting capacity to cover all internal moving parts
- Heat resistant and does not evaporate
- 500ml

PEGA-100111

£11.10



T1 Oil Spray

- Universal cleaning and lubricating oil for the preservation and disinfection of dental instruments
- 250ml

Dentsply Sirona



SIR-5901665

£39

F1 Service Oil Spray MD

- Superior lubrication
- Fully sterilizable, absolutely resistant to hydrolysis
- F1 Oil works with turbines, straight and contra-angle handpieces, surgical instruments, air scalers and air motors
- 400ml

WH-10940021

£36.95



Oil Spray Nozzle - Fixed Connection

- Spray Cap With Nozzle
- Chuck system for turbines and contra-angle handpieces as well as turbine handpieces and air motors with fixed connection

WH-02036100

£9



Oil Spray Nozzle - Slow Handpieces

- Spray cap for W&H F-1 service oil
- For slow handpieces

WH-02038200

£12



Oil Spray Nozzle - Fast Rotoquick Connection

- Spray cap with nozzle for all W&H turbines, air motors and scalers with ROTO QUICK connection

WH-02229200

£20



Pana Spray Oil

- PANA SPRAY Plus is a specially formulated, highly effective general handpiece and air motor lubricant.

NSK



NSK-Z182100

£29.43

iCare Maintenance Oil

- Specially formulated maintenance oil provides for the highest levels of cleaning & lubricating of handpieces
- 1000ml

NSK



NSK-Z016130

£84.39

Woodpecker O - Star

The Ultimate Curing Expert

- Performs light curing in 1 second
- Longer wavelength: 385 - 515 nm
- It includes a light intensity meter which is located in the charging base
- It has 7 different work modes: normal, high intensity, turbo, ortho, soft, pulse and verification
- Light intensity: up to 3000mw/cm²
- Small, rotating head allows you to reach all angles of the mouth and reach the back teeth without difficulty

WOOD-OSTAR

£399



WOODPECKER

Woodpecker iLED

Curing Light

- Ergonomic curing light with 360° rotating head
- Fast - Cures up to 2mm in just one second
- Light intensity up to 2300mw/cm² with constant power output
- Power input: AC100V-240V
- Light output intensity: 1000 mW/cm²~2500mW/cm²
- Time Sets: 5-20 Seconds
- Operation Mode: turbo mode, normal mode

WOOD-ILED

£199



WOODPECKER

Woodpecker LED.H ORTHO

Specially designed for Orthodontics

- Specially designed for orthodontics. In just 3 seconds the brace solidifies
- Short, powerful cycles without overheating
- Constant work without losing intensity, for a complete curing of the resins
- Light and easy to use, with a rechargeable battery
- The new bevelled light guide design can easily reach any tooth and radiate vertically at any angle
- Includes a light intensity meter

WOOD-LED.H.ORTHO

£195



WOODPECKER

Woodpecker LED.H

Specially designed for Orthodontics

- Woodpecker LED-H Light Curing Lamp provides reliable light curing of any type of material
- It is a light weight lamp, which will allow you to use it during a large part of your professional daily activity, without causing fatigue
- Three working modes: Normal, ramp or intermittent
- High-capacity battery with charge duration for 200 continuous light-curing cycles

WOOD-LED.H

£160



WOODPECKER

Sine Wave Reservoir Scaler

- Ultrasonic cleaning device with a handpiece type EMS without light
- Dismountable and sterilizable table-top model
- Modern design with touch control and waterproof
- Scaler, perio and endo function
- Modulates sinus waves, allowing the best treatment with maximum comfort

Includes:

- 7 tips: 1 x G4 supragingival scaler tip
- 1 x G12 supragingival scaler tip
- 1 x P1 Periodontal scaler tip
- 1 x P3 Periodontal scaler tip
- 1 x P8 Periodontal scaler tip
- 2 x E6 Endodontic scaler tip

WOOD-U600

£595



WOODPECKER

Ultrasonic UDS-N2 Scaler

- Built in ultrasonic scaler (compatible with A-dec and Belmont units)
- The handpiece is detachable and can be autoclaved to 135 °
- The UDS N2 instrument is compatible with the EMS system compatible tip and inserts

WOODPECKER

£450

WOOD-UDS-N2



DTE-D7 Piezo Scaler

Satelec Type

- Assembled with excellent piezo-ceramic, the handpiece has high and lasting power and generates little heat.
- Eliminates dental calculus, bacterial plaque and tea scales on the teeth quickly without hurting the gum or the enamel.
- Light weight, easy to operate.
- Compatible with Satelec hand piece and tips.



£399
DTE-D7

Replacement Handpieces

- Woodpecker/EMS Type WOOD-HW-3H £145
- Woodpecker U600 Handpiece WOOD-HW-6 £145
- Woodpecker/EMS Type LED WOOD-HW-5L £160
- DTE/Acteon Satelec Type LED DTE-HD-7H £145

Replacement Wrenches

- Woodpecker/EMS Type WOOD-TW-IL £10.50
- Sirona Type WOOD-TS-IL £13
- DTE/Acteon Satelec Type DTE-TD-IL £10.50
- Satelec Type DTE-TD-2L £13

Replacement Tips

- | | |
|--------------------|---------------------|
| Woodpecker/EMS | DTE/Acteon/Satelec |
| Perio P1 WOOD-P1 | Perio PD1 DTE-PD1 |
| Perio P3 WOOD-P3 | Perio PD3 DTE-PD3 |
| Perio P4 WOOD-P4 | Scaling GD1 DTE-GD1 |
| Scaling G1 WOOD-G1 | Scaling GD4 DTE-GD4 |
- £18**