

I ♥ KaVo DEALS

Black +
Friday

Special

Extended offer - get it now!



The limited I ♥ KaVo DEALS with prices to fall in love with,
once a month for two popular KaVo products.

EXPERTmatic LUX E25 L & EXPERTtorque LUX E680 L.

For this week only, take advantage of phenomenal prices for the

KaVo EXPERTmatic LUX E25 L (KAVO-E25L) for a **special price** of the I ♥ KaVo promotion of **€ 561 each + VAT**,

Normally €906 + VAT

KaVo EXPERTtorque LUX E680 L (KAVO-E680L) for a **special price** of the I ♥ KaVo promotion of **€ 442 each + VAT**,

Normally €764 + VAT


This price is valid from **20.-28. November 2023**.





The KaVo EXPERT Series: Rely on KaVo's advanced performance.

- Small head size for easy, efficient access to the preparation area
- High-tech FG chuck for optimum longevity
- 4 or 3-port spray with ultra-fine spray wetting for optimum cooling
- Push-button chuck system with a retention power of up to 30N for secure retention of the bur
- Glass rod light conductor with up to 25,000 lux for reliable detection of the preparation point
- Ceramic ball bearings for excellent quiet and low-vibration operation
- Low volume to eliminate unpleasant noises for patients and treatment personnel
- Optimum access thanks to the ergonomic 100°/19° head/knee angle combination of the EXPERT turbine
- Thermo-disinfectable and sterilisable up to 135 °C
- 24-month warranty



KaVo EXPERTmatic LUX E25 L

- Fast, precise work thanks to small head
- Exclusive ceramic ball bearings
- Internally channelled spray water and spray air
- CoolHead technology

KAVO-E25L



KaVo EXPERTtorque LUX E680 L

- Head height of 13.6 mm
- Maximum power: up to 29 watts
- Low volume: average 59 dB(A)
- With glass rod light conductor

KAVO-E680L



KAVO
Dental Excellence