

I ♥ KaVo DEALS

Black +
Friday

Special

Extended offer - get it now!



The limited I ♥ KaVo DEALS with prices to fall in love with,
once a month for two popular KaVo products.

EXPERTmatic LUX E25 L & EXPERTtorque LUX E680 L.

For this week only, take advantage of phenomenal prices for the

KaVo EXPERTmatic LUX E25 L (KAVO-E25L) for a **special price** of the I ♥ KaVo promotion of **£474 each + VAT**,
Normally £788 + VAT

KaVo EXPERTtorque LUX E680 L (KAVO-E680L) for a **special price** of the I ♥ KaVo promotion of **£373 each + VAT**,
Normally £665 + VAT


This price is valid from **20.-28. November 2023**.





The KaVo EXPERT Series: Rely on KaVo's advanced performance.

- Small head size for easy, efficient access to the preparation area
- High-tech FG chuck for optimum longevity
- 4 or 3-port spray with ultra-fine spray wetting for optimum cooling
- Push-button chuck system with a retention power of up to 30N for secure retention of the bur
- Glass rod light conductor with up to 25,000 lux for reliable detection of the preparation point
- Ceramic ball bearings for excellent quiet and low-vibration operation
- Low volume to eliminate unpleasant noises for patients and treatment personnel
- Optimum access thanks to the ergonomic 100°/19° head/knee angle combination of the EXPERT turbine
- Thermo-disinfectable and sterilisable up to 135 °C
- 24-month warranty



KaVo EXPERTmatic LUX E25 L

- Fast, precise work thanks to small head
- Exclusive ceramic ball bearings
- Internally channelled spray water and spray air
- CoolHead technology

KAVO-E25L



KaVo EXPERTtorque LUX E680 L

- Head height of 13.6 mm
- Maximum power: up to 29 watts
- Low volume: average 59 dB(A)
- With glass rod light conductor

KAVO-E680L

Promo prices apply exclusively to individual instruments, not to bundles and cannot be combined with other promotions. The promo prices quoted do not include VAT and refer to the list prices of KaVo Dental to dealers. List prices of the dealer remain unaffected. Dealers are free to decide whether and in what form they participate in the promo. The November Deals are valid from 20.-28.11.2023.



KAVO
Dental Excellence