

## Surf the canal with confidence







## Keeping the WAVEONE® simplicity, with a golden host of benefits




## Reinforces patient safety

- Primary WAVEONE® GOLD file is 50% more resistant to cyclic fatigue than WAVEONE® Primary file
- Reduced screwing effect compared to standard rotary systems

# Cyclic fatigue resistance 50% WAVEONE® NIT i systems Primary GOLD Primary GOLD Primary


## Covers a wider range of canal morphologies

- Enhanced file flexibility thanks to the GOLD technology
- Extended size range (Small, Primary, Medium, Large)



## Shortens the shaping time

- Takes the cutting efficiency to a higher level
- A single file per treatment also translates into time savings for canal shaping and irrigation



## Engineered with the latest

## reciprocating technology

DENTSPLY MAILLEFER has teamed up with the best endodontic professionals to develop wave.one

#### The design team speaks

Dr C. Ruddle, Dr J. Webber, Dr W. Pertot, Dr S. Kuttler Shaping canals has never been easier.

Since its launch in 2011, Dentsply Maillefer's WAVEONE® has set a new standard in root canal shaping. All around the world general dentists and endodontists have voiced their appreciation of the simplicity WAVEONE® brings to their efforts to provide quality endodontic treatment to their patients.


We have been working on a new and exciting WAVEONE® project to bring you WAVEONE® GOLD.

WAVEONE® GOLD brings to the endodontic procedure a metallurgically advanced single-file technique for shaping canals. Further, we optimized the tip diameters, tapers, and the cross-section to produce a file that really improves safety, efficiency and flexibility when preparing canals. WAVEONE® GOLD takes simplicity and treatment success to another level.










## Simplicity that enables treatment focus

## The majority of root canals can be shaped with a single wave one Primary file





### Ready-to-use + Single use








Most dental practitioners have the twofold aim of reducing overall treatment costs while improving safety and efficiency. With this in mind, WAVEONE® GOLD files are supplied in sterilized blister packs. Dentists will benefit from:

- Simplicity, thanks to the ready-to-use sterile packaging
- Increased safety margin, by always using files with maximum fatigue resistance potential
- Reduced risk of cross-contamination, which can occur when instruments are reused
- Maximized cutting efficiency

## One complete solution

With wave • one , Dentsply Maillefer promotes comprehensive treatment solution and optimizes predictably successful outcomes



#### Files



#### small reciprocating file

3x small (020 07) - 21mm - A 0753 221 GOS 03

**3x** small (020 07) - 25mm - A 0753 225 GOS 03

**3x** small (020 07) - 31mm - A 0753 231 GOS 03

6x small (020 07) - 21mm - A 0756 221 GOS 03

**6x** small (020 07) - 25mm - A 0756 225 GOS 03 **6x** small (020 07) - 31mm - A 0756 231 GOS 03

primary reciprocating file

3x primary (025 07) - 21mm - A 0753 221 GOP 03

**3x** primary (025 07) – 25mm – A 0753 225 GOP 03

**3x** primary (025 07) – 31mm – A 0753 231 GOP 03

6x primary (025 07) - 21mm - A 0756 221 GOP 03

**6x** primary (025 07) - 25mm - A 0756 225 GOP 03

**6x** primary (025 07) - 31mm - A 0756 231 GOP 03

#### medium reciprocating file

**3x** medium (035 06) - 21mm - A 0753 221 GOM 03

**3x** medium (035 06) - 25mm - A 0753 225 G0M 03

**3x** medium (035 06) - 31mm - A 0753 231 GOM 03

**6x** medium (035 06) - 21mm - A 0756 221 GOM 03

6x medium (035 06) - 25mm - A 0756 225 GOM 03

**6x** medium (035 06) - 31mm - A 0756 231 GOM 03

#### large reciprocating file

**3x** large (045 05) - 21mm - A 0753 221 GOL 03

3x large (045 05) - 25mm - A 0753 225 GOL 03

**3x** large (045 05) - 31mm - A 0753 231 GOL 03

6x large (045 05) - 21mm - A 0756 221 GOL 03

**6x** large (045 05) - 25mm - A 0756 225 GOL 03 6x large (045 05) - 31mm - A 0756 231 GOL 03

Oven

THERMAPREP® 2

A 0176 115 002 00 - (115V)

A 0176 230 002 00 - EU (230V)

A 0176 230 002 33 - UK (230V)

#### assortment reciprocating file

1x small/primary/medium/large - 21mm

A 0754 221 G9A 03

**1x** small/primary/medium/large – 25mm

A 0754 225 G9A 03

1x small/primary/medium/large - 31mm

A 0754 231 G9A 03

## Paper points

#### Box of 180x sterile paper points

small - A 175W 000 00S 03

primary - A 175W 000 00P 03

medium - A 175W 000 00M 03

large - A 175W 000 00L 03

## Gutta percha points

#### Box of 60x variable tapered gutta percha points

small - A 175X 000 00S 00

primary - A 175X 000 00P 00 medium - A 175X 000 00M 00

large - A 175X 000 00L 00

#### Blister of 6x obturators

small obturators - A 1750 B00 00S 00 primary obturators - A 1750 B00 00P 00

medium obturators - A 1750 B00 00M 00

large obturators - A 1750 B00 00L 00

#### Box of 30x obturators (5 blisters)

small obturators - A 1750 POO OOS OO

primary obturators - A 1750 P00 00P 00

medium obturators - A 1750 POO 00M 00

large obturators - A 1750 POO OOL OO

#### Kits

#### A10541 X-SMART IQ™ Starter Kit:

1 Motor with contra-angle

1 Pack of 15 Primary WAVEONE GOLD

1 Pack of 16 Assorted WAVEONE GOLD

1 Blister of 3 PROGLIDER

1 iPad Protective Case

#### A10358 X-SMART® Plus Starter Kit:

1 Motor with contra-angle

1 Pack of 15 Primary WAVEONE GOLD

1 Pack of 16 Assorted WAVEONE GOLD

1 Blister of 3 PROGLIDER



Motors

X-SMART IQ™ Motor - A1054

X-SMART® Plus Motor - A1032



